

Sales of this product have ended.

Terahertz Wave Spectroscopy and Analysis Platform

Full Coverage of Applications From R&D to Industrial Testing




Optimal for a wide range of terahertz research and applications

The TAS7500TS Terahertz Wave Spectroscopy and Analysis Platform consists of an optical fiber laser module and a data acquisition unit. Flexible source selection and source/detector placement allow the user to build customized configurations




Basic Configuration



▶ TAS7500TS/TAS7400TS : Superior Performance


Excellent Reproducibility

Advantest's proprietary low-jitter fiber laser and analog signal analysis technology enable superior refractive index (phase) measurement repeatability of within $\pm 0.03\%$.



Broadband Terahertz Spectroscopy

In addition to standard type TAS1110, two specialized source modules - the low frequency type TAS1120 and the broadband type TAS1130 - expand the bandwidth coverage of the platform to serve a broad array of applications.



Application

Spectroscopic analysis of terahertz optical devices and chemical materials Materials research and development in the sub-terahertz range

Spectroscopy and analysis in the terahertz region up to 7 THz

Target Materials Polymers/polymorphic crystalline forms/biomaterials/ controlled substances/ metamaterials/food/crops/ construction materials

Dielectric materials/paint/ plastics/construction materials Aqueous solutions/liquid solutions/oxides/biomaterials/ functional devices/graphene/ resinous materials

Frequency Range

 $0.1 \sim 4 \, \text{THz}$

 $0.03 \sim 2$ THz

 $0.5 \sim 7$ THz

THz source Module

TAS1110

TAS1120

TAS1130

THz Detector Module **TAS1230**

Terahertz Optical Sampling System

TAS7500TS 1ms/scan

TAS7400TS 200 ms/scan

FFATURES

- The system uses two channels of ultra short pulse lasers (1550 nm) with either biased output (for THz generation) or signal input (for THz detection). Advantest's unique optical sampling method, utilizing phase-modulated dual-laser-synchronized control technology without a mechanical optical delay line, enables extremely high speed terahertz spectroscopy.
- The basic configuration has a single channel input and output, and an optional second channel with software controlled duplexer can
- Spectroscopic functionality is implemented by connecting optional fiber-coupled THz source and detector modules.
- Ethernet remote control option enables support for remote programming of THz measurement and analysis functions and peripheral equipment.



Terahertz Source Modules

► TAS1110 → TAS1120 → TAS1130


Standard: 0.1~4 THz

Low frequency type: $0.03\sim$ 2 THz

Broadband type: $0.5\sim7$ THz

FEATURES

- Easy to use all-in-one THz source modules complete with a fiber pigtailed compact housing.
- ▼ TAS1110 and TAS1120 consist of a photoconductive antenna and a hyper-hemispherical silicon lens. Current monitor and thermistor on the built-in bias circuit deliver stable output
- intensity regardless of environmental temperature changes. TAS1130 is a Cherenkov THz source module that utilizes non-linear optical crystal LiNbO3 waveguides.



•THz frequency comb generation



Terahertz Detector Module

▶ TAS1230

FEATURES

- Easy to use all-in-one photoconductive THz detector complete with hyper-hemispherical silicon lens and fiber
- pigtailed compact housing. THz waves are sampled by introducing 1550 nm ultra short
- laser pulses. Built-in trans-impedance (current-voltage transfer) amplifier with 500 kHz bandwidth simplifies implementation.






Usage: •THz wave generation for time domain spectroscopy (TDS), etc.

System Configurations with the Optical Bench Unit Option

Supports Measurement & Imaging via Transmittance and Reflectance Methods


- Can measure and render 2D images of specimens up to a maximum size of 150mm x 150mm
- Remote programming option supports Visual Basic, Excel VBA, C#, and LabVIEW programming languages. Calculation libraries enable users to easily create programs including analysis of acquired data.
- Sample programs are provided for each language, enabling users to build system configurations with linked peripheral equipment.



Dual Channel Monitoring Option

Supports highly stable measurement for long measurement times.

Dual beam path architecture enables highly stable measurement over long periods. By using this option, the system can monitor the amplitude and phase fluctuation during the measurement using the second detector and compensates automatically.



Specifications

■ TAS7500TS/TAS7400TS Terahertz Optical Sampling System

Items		Specification						
			TAS7500TS	TAS7400TS				
System configuration		Dual ultra-short pulsed optical fiber laser (for THz generation and detection) measurement unit, analysis unit, controller & analysis PC						
Excitation laser	Center wavelength	1550 nm						
	Output power	≥ 20 mW (Option: Max. ≥50 mW)*1						
	Pulse width	≤ 50 fs (Using 1.5 m fiber)*1						
	Repetition rate	50 MHz \pm 200 Hz						
	Optical output port	THz generation: 1 port, THz detection: 1 port (can be optionally increased to 2 of each)						
Measurement specifications	Measurement method	Terahertz optical sampling method (phase modulation method)						
	Time resolution	2 fs						
	Frequency resolution	3.8 GHz	7.6 GHz	61.0 GHz	1.9 GHz	7.6 GHz		
	Scan range	262 ps	131 ps	16 ps	524 ps	131 ps		
	Throughput	16 ms/scan	8 ms/scan	1 ms/scan	200 ms/scan			
	Frequency accuracy	±10 GHz*1						
General specifications	Usage environment	Temperature range: +10°C to +30°C Relative humidity: 80% or less (no condensation)						
	Storage environment	Temperature range: - 10° C to + 50° C Relative humidity: 80% or less (no condensation)						
	Power	AC100V(100-120)/200V(220-240)±10%, 50/60 Hz, 250VA (not including analysis PC)						
	Size/Weight	Analysis unit: 430 (W) \times 540 (D) \times 330 (H) mm / Analysis unit: 30 kg or less Measurement unit: 430 (W) \times 240 (D) \times 220 (H) mm / Measurement unit: 14 kg or less						

^{*1:} At temperatures of:23°C±5°C

■ TAS1110/1120/1130 Terahertz Source Module

Items			Notes			
		TAS1110 (Standard)	TAS1120 (Low frequency Type)	TAS1130 (Broadband type)		
Generation method		Photoconductive antenna		Cherenkov type THz wave generation using LiNbO ₃ waveguide		
Bandwidth (SNR=1)*2*3		0.1 to 4 THz	0.03 to 2 THz	0.5 to 7 THz	Input power: 20 mW (TAS1110 / TAS112	
Dynamic range*2*4 (Peak level)	with the TAS7500TS*5	≥ 70 dB	≧ 60 dB	≥ 70 dB	50 mW (TAS1130) Laser pulse width: 50 fs	
	with the TAS7400TS*6	≥ 60 dB	≥ 50 dB	≥ 60 dB		
Input optical fiber connector		ø 3 mm 1550 nm Polarization maintaining fiber			Length: 1.5 m	
Size		55 mm × 20 mm × 20 mm		43 mm × 24 mm × 21 mm	Without fiber pigtail	

■ TAS1230 Terahertz Detector Module

Items		Specification	Notes	
Detection method		Photoconductive antenna		
Bandwidth (SNR=1)		0.1 to 4 THz (Generator: TAS1110) 0.5 to 7 THz (Generator: TAS1130)		
Dynamic range (Peak level)	with the TAS7500TS	≥ 70 dB (Detector: TAS1110, TAS1130) ≥ 60 dB (Detector: TAS1120)	Input power:20 mW Laser pulse width: 50 fs	
	with the TAS7400TS	≧ 60 dB (Detector: TAS1110, TAS1130) ≧ 50 dB (Detector: TAS1120)		
TIA sensitivity		$9.7 \times 10^6 \text{ V/A}$		
TIA bandwidth		500 kHz	-3 dB bandwidth	
Input optical fiber connector		ø 3 mm 1550 nm Polarization maintaining fiber	Length: 1.5 m	
Size		55 mm × 20 mm × 20 mm	Without fiber pigtail	

ADVANTEST

https://www.advantest.com/

ADVANTEST CORPORATION

• Futher Inquiries

E:mail: info_t@advantest.com Phone: +81-480-72-6300

Advantest Corporation

1-5, Shin-tone, Kazo-shi, Saitama 349-1158 Japan

^{*2:} Detector: TAS1230 *3: In case of measuring sample with the TAS7500TS, at frequency resolution: 7.6 GHz.

*4: The peak level frequency varies in each system, and the dynamic range on each frequency varies in each system. *5: Number of integration: 262144 *6: Number of integration: 16384